**[ENGLISH] INSTRUCTIONS FOR USE****I. INTRODUCTION**

KATANA AVENCIA Block is CAD/CAM blocks of composite resin material that is milited to fabricate dental restorations, including inlays and crowns.

II. INDICATIONS

KATANA AVENCIA Block is indicated for the following uses:
Fabrication of inlays, onlays, veneers and full contour crown restorations using a dental CAD/CAM system.

III. CONTRAINDICATIONS

- Patients with a history of hypersensitivity to methacrylate monomers
- This product should NOT be used when malocclusion, clenching or bruxism conditions are present.

IV. PRECAUTIONS

- Safety precautions**
1. This product contains substances that may cause allergic reactions. Avoid use if the product is patients with known allergies to methacrylate monomers or any other components.
2. The patient demonstrates a hypersensitivity reaction, such as rash, eczema, features of inflammation, ulcer, swelling, itching or numbness, discontinue use of the product and seek medical attention.
3. Avoid direct contact with the skin and/or soft tissue to prevent hypersensitivity. Wear gloves or take appropriate precautions when using the product.
4. Single use only. Do not reuse the product.
5. Exercise caution to prevent the patient from accidentally swallowing the product.
6. If the instruments of this product are damaged, protect yourself from any danger and immediately discontinue their use.
7. Do not use the product if there is a crack noticed after removing it from the package. Do not use if you find a crack in the final prosthesis.
8. Dispose of this product as a medical waste to prevent infection.

2. Handling and manipulation precautions

- Precautions when using the product.
1. The product must not be used for any purposes other than specified in [INDICATIONS].
2. The use of this product is restricted to dental professionals.
3. Use safety equipment, such as a local vacuum unit, protective shields, and an appropriately approved dustproof mask, to prevent the patient from inhaling dust when shaping, finishing or polishing the product.
4. For bonding, use a dental cement, e.g. PANAVIA VS or PANAVIA SA Cement (Universal). Before bonding, sandblast the adherent surface and then condition it with an etching agent (e.g. K-ETCHANT Syringe) and with a silane coupling agent (e.g. CLEARFIL CERAMIC PRIMER PLUS), according to the Instructions for Use.
5. Before applying additional resin cement, characterizing or glazing the surface, treat the adherent surface by etching it with an etching agent and then condition it with a silane coupling agent, according to the Instructions for Use.
6. Shape the outlines of inlays or onlays with appropriate smooth curves and locate the margins so that there is no contact with the opposing teeth.
7. The margins of crowns should be prepared with a deep chamfer.
8. When preparing the abutment of the crown, avoid the following:
• Deep shoulder or "J" margin
• Knife Edge
• Rough Margin
• Parallel Axis
• Undercut
• Groove
• Sharp Incisal Tip
9. When using this product to fabricate a restoration, be sure to attain the thicknesses listed in the following table:

Restoration		Thickness
Anterior Crown	Incisal	>1.5 mm
	Circumferential wall	1.0-1.2 mm
	Margin	>0.8 mm
Posterior Crown	Fissure / Cusps / Circumferential wall	>1.5 mm
	Margin	>1.0 mm
	Fissure	>1.0 mm
Inlay	Isthmus	>1.5 mm
	Fissure	>1.0 mm
	Cusps / Isthmus	>1.5 mm
Onlay	Incisal	0.5-1.5 mm
	Labial wall	>0.6 mm
	Margin	>0.4 mm

- This product is a CAD/CAM composite resin material. When milling the product, use the appropriate CAD/CAM milling machine and follow the Instructions for Use.
- Do not bake this product in a porcelain furnace or similar item.

3. Storage precautions

- The product must be used by the expiration date indicated on the package.
2. The product must be stored at 2-25°C/36-77°F when not in use. Store in a cool and dry place. Keep away from direct sunlight.
3. The product must be stored in a proper place where only dental practitioners can access.

V. COMPONENTS**1. Shades**

KATANA AVENCIA Block is available in the following 3 shades:
A2 LT, A3 LT, A3.5 LT

2. Sizes

KATANA AVENCIA Block is available in the following 2 sizes:
Size 12 (10x12x15 mm)
Size 14L (14.5x14.5x18 mm)

3. Metallic holder

KATANA AVENCIA Block is 2 kinds of metallic holder available: Universal and CEREC.

4. Components

Please see the outside of the package for contents and quantity.

5. Ingredients

- Principal ingredients:
- Mixed filler with colloidal silica and aluminum oxide
 - Cured resins consisting of methacrylate monomer
 - (Copolymer of Urethane dimethacrylate and other methacrylate monomers)
 - Pigments

VI. CLINICAL PROCEDURES**1. Preparation**

Prepare the restoration area in the usual manner.

2. Designing and processing the prosthesis
Select the product shade and size that suits the case. Design and process it using a CAD/CAM system. See the Instructions for Use.

3. Shaping and polishing
Shape the prosthesis in the usual manner and polish it with silicon points and other available methods.

4. Sandblasting
Roughen the adherent surface by blasting with aluminum powder (30-50 µm) using air pressure of 0.1-0.2 MPa (1-2 kgf/cm², 14-29 psi).

5. Bonding

After a final fit, treat the adherent surface with an etching agent and a silane coupling agent according to the Instructions for Use. Use a resin cement to bond the prosthesis.

[WARRANTY]

Kuraray Noritake Dental Inc. will replace any product that is proven to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

[NOTE]

If a serious accident attributable to this product occurs, report it to the manufacturer's authorized representative shown below and the regulatory authorities of the country in which the user/patient resides.

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"KATANA" is a trademark of NORITAKE CO., LIMITED.
"AVENCIA", "PANAVIA" and "CLEARFIL" are trademarks of KURARAY CO., LTD.
A "KATANA" a NORITAKE CO., LTD. védjegye.
A "AVENCIA", a "PANAVIA" és a "CLEARFIL" a KURARAY CO., LTD. védjegyei.
Kuraray Europe GmbH (Importer)