

Eztetic™ 3.1 mmD Implant System

Product Catalog And Surgical Technique



How To Order

How To Order

To order, call Zimmer Biomet Dental Customer Service between 7:00 a.m. and 7:00 p.m. EST, Monday through Friday, or fax your order:

U.S. Toll-Free: 1 (800) 342-5454
Tel: +1 (561) 776-6700
Fax: +1 (561) 776-1272

Overnight delivery is available, if requested. Overnight delivery orders must be placed by 7:00 p.m. EST.

Specifications and prices are subject to change without notification.

Please note that not all products are registered or available in every country/region. Please check with a Zimmer Biomet Dental representative for availability and additional information.

Images within this catalog not to scale.

Zimmer Biomet Dental accepts VISA, MasterCard and American Express.

Payment terms are credit card or net 30. All payments must be in U.S. Dollars (no foreign currency accepted).

Return Policy

A purchaser may return any Zimmer Biomet Dental product within 180 days of invoice date. Zimmer Biomet Dental products returned within 90 days of the invoice date will be accepted without any restocking fee to the purchaser. Zimmer Biomet Dental products returned between 91 and 180 days of the invoice date will be subject to a 15% restocking fee.

Exchange Policy

Zimmer Biomet Dental products may be exchanged within 365 days of the invoice date, without a restocking fee, for other Zimmer Biomet Dental products within the same product family as determined by Zimmer Biomet Dental (e.g., implant for implant, abutment for abutment, membrane for membrane, etc.), provided that if the replacement products are of greater value than the returned products, the purchaser will be charged the price differential.

The value of a replacement product will be determined based on the following in order of priority: promotional discount applied to purchase the returned item so long as the original promo code or pricing contract also covered the replacement product; or list price. Freight charges will apply to exchanged product shipments.

Conditions Applicable to All Returns & Exchanges

Any returned or exchanged Zimmer Biomet Dental product must be (1) returned in its original, unopened package (including autoclavable product); (2) a non-expired and currently offered Zimmer Biomet Dental product (expired and discontinued products are excluded); (3) returned freight prepaid to 4555 Riverside Drive, Palm Beach Gardens, Florida 33410


(Attention: Returns Department) accompanied by the Zimmer Biomet Dental Return Authorization Number provided by Zimmer Biomet Dental's Customer Service Department and (4) received by the Returns Department within the time frames specified above. Zimmer Biomet Dental products not meeting these criteria will not be accepted for return or exchange. Returned and exchanged product lot numbers will be checked to ensure compliance with the time limitations referenced in this policy. Zimmer Biomet Dental reserves the right, in its sole discretion, to (a) decline any return or exchange if product with alternative lot numbers is provided in place of product purchased within the time limits herein, or (b) terminate the eligibility of any customer from this courtesy return policy in the event of suspected abuse or fraud.



Table Of Contents

Overview	1
• General Information	1

Implant Design And Specifications	3
• Implant Platform	3
• Implant Placement Guidelines	4
• Indications For Use	4



Eztetic Implants	5
• Implants	5
• Healing Collars	5



Eztetic Prosthetics (NP Platform)	6
• Impression Transfer Components	6
• Provisional Restoration Components	6
• Cement-Retained Components	6
• “Cast To” Gold Abutment Components	7
• Overdenture Attachments	7
• Prosthetic Accessories And Tools	8



Eztetic Instrumentation	9
• Color Coding	9
• NP Surgical Module	9
• Replacement Instrumentation	9
• Surgical Kit	9



Surgical Technique	12
• Drilling Sequence Guidelines	12
• Surgical Procedure	13




References	14
-------------------	----

Overview

1


General Information

The 3.1 mmD Eztetic Implant offers a narrow, powerful solution for demanding anterior spaces. This innovative implant is designed to provide initial stability¹, allowing immediate restorations when clinically appropriate. The conical, Double Friction-Fit™ Connection with platform switch, combined with narrow platform (NP) prosthetics, are designed for crestal bone maintenance² and optimum aesthetics by accommodating maximum soft-tissue volume. The NP Surgical Module conveniently snaps into the Tapered Screw-Vent® Surgical Kit for streamlined surgical procedures.



Restorative Profile For Aesthetic Emergence:

Implant-abutment connection along with a contour abutment profile are designed to provide space for soft-tissue and aesthetic emergence of the restoration





Strength³ For Long-Lasting Aesthetics

Implant design and a conical Double Friction-Fit Connection are combined for exceptional strength, reduced micromovement and microleakage

Primary Stability¹ For Immediate Aesthetics

Tapered implant geometry combined with dedicated soft- and dense-bone surgical protocols are designed for high primary stability in all types of bone



MTX[®] Microtexture
At 2000x Magnification

MTX Surface For Increased Bone Apposition^{3,4}

The MTX Microtextured Surface has been documented to achieve high levels of bone-to-implant contact and successful clinical results under conditions of immediate loading


Coronal Options For Bone Level Maintenance

The coronal microgrooves are designed to preserve crestal bone.⁵

Two coronal surface configurations are available:

- Full MTX Microtexturing with MTX Crestal Microgrooves (Model CT)
- 0.5 mm Machined Collar with MTX Crestal Microgrooves (Model CM)

Eztetic Dental Implant - Fully Textured With Microgrooves And 0.5 mm Machined Collar With Microgrooves



The 3.1 mmD Eztetic Dental Implants have a MTX microtextured or 0.5 mm machined coronal aspect, followed by 1.8 mm of the MTX Surface with microgrooves. The six microgrooves are circumferential with a depth of 0.06 mm and peak-to-peak width of 0.3 mm. Triple-lead threads begin immediately after the microgrooves and continue to the apex. The degree of body taper varies between 1° and 4°, depending on implant length.

Conical, Double Friction-Fit Technology

The 3.1 mmD Eztetic Implants feature a conical Double Friction-Fit Connection designed to reduce micromovement and microleakage via a precise implant-abutment interface, aiding in crestal bone maintenance.

Figure 1A The implant-abutment connection offset along with a Contour Abutment are designed to provide space for soft-tissue and aesthetic emergence of the restoration.






Figure 1B The 3.1 mmD Eztetic Implants feature a conical Double Friction-Fit Connection designed to reduce micromovement and microleakage via a precise implant-abutment interface, aiding in crestal bone maintenance.




Implant Platform

The implant platform diameter is measured across the most coronal part of the implant. The 3.1 mmD Eztetic Implant features a 2.9 mmD prosthetic platform. A 1.3 mm deep, 17° internal cone extends from the outermost diameter (2.9 mmD) of the implant platform to the internal hex of the implant. The internal hex is 2.1 mm flat-to-flat with a depth of 1.7 mm (Fig. 2A and 2B). The 3.0 mm deep conical connection is designed to distribute stresses deep into the implant and away from the crestal bone to aid in crestal bone maintenance.



Implant Placement Guidelines

For Pre-Surgical Planning Instructions and Anatomical Criteria guidelines please refer to the Tapered Screw-Vent Surgical Manual.



Collar Depth Adjustment

The 3.1 mmD Eztetic Dental Implants are designed to be placed at bone level or slightly below bone level, according to clinical preference and patient anatomy.

Indications For Use

The 3.1 mmD Eztetic Dental Implants are designed for use in the anterior maxilla or mandible for immediate loading or for loading after a conventional healing period. Implants may be used to replace one or more missing teeth. Immediate loading is indicated when there is good primary stability and an appropriate occlusal load.

The 3.1 mmD Eztetic Dental Implants may be placed immediately following an extraction or loss of natural teeth provided there is sufficient volume of alveolar bone to minimally support the implant (minimum 1.0 mm circumferential and 2.0 mm apical).

The 3.1 mmD Eztetic Dental Implants should be splinted to additional implants when used in the pre-molar region and should not be used in the molar region. For more information, please review the full instructions for use.

Implant Platform Color-Coding

3.1 mmD Eztetic Implants are available with a 2.9 mmD prosthetic platform.

Implant Diameter	Implant Platform	Color-Coding
3.1 mmD	NP (2.9 mmD)	● Light Blue

NP = Narrow Platform

Abutment Emergence Profile Compatibility

Abutment Emergence Profile*	Color-Coding
4.5 mmD	● Tan

*For Contour components.

NOTE: Contour Abutments with 3.7 mmD Emergence Profile are not compatible with 3.5 mmD Contour Restorative Components.

Eztetic Dental Implants, MTX Surface, Fully Textured With Microgrooves

Includes Healing Screw (CCSNP)



Catalog Numbers						
Implant Diameter	Implant Platform	8.0 mmL	10 mmL	11.5 mmL	13 mmL	16 mmL
3.1 mmD	● NP (2.9 mmD)	CT318	CT3110	CT3111	CT3113	CT3116

Eztetic Dental Implants, MTX Surface, 0.5 mm Machined Collar With Microgrooves


Includes Healing Screw (CCSNP)



Catalog Numbers						
Implant Diameter	Implant Platform	8.0 mmL	10 mmL	11.5 mmL	13 mmL	16 mmL
3.1 mmD	● NP (2.9 mmD)	CM318	CM3110	CM3111	CM3113	CM3116

Implant Fatigue Strength¹


The 3.1 mmD Eztetic Implants achieved 43% higher fatigue strength compared to selected competitive implants of similar diameters.¹



All Products were tested in increments of 5.

Insertion Torque¹


The 3.1 mmD Eztetic Implants achieved high insertion torque.¹



Benchtop engineering test utilizing a dense bone substrate.¹

Torsional Yield Strength¹

The 3.1 mmD Eztetic Implant interface withstood higher torsional forces than the selected competitors.¹



Benchtop engineering test utilizing the implants and their corresponding drivers.¹

Surgical Healing Screw, Replacement



Catalog Numbers	
Implant Platform	Catalog No.
● NP (2.9 mmD)	CCSNP

Healing Collars



Catalog Numbers				
Implant Platform	Emergence Profile	Cuff Height		
		1.5 mm	3.0 mm	4.5 mm
● NP (2.9 mmD)	3.7 mmD	CHCNP31	CHCNP33	CHCNP34
● NP (2.9 mmD)	4.5 mmD	CHCNP41	CHCNP43	CHCNP44

Indirect Transfers (Closed-Tray Procedure) Color-coded by implant platform.

Includes a Retaining Screw (CASLT)



Catalog Numbers				
Implant Platform	Emergence Profile	Cuff Height		
		1.5 mm	3.0 mm	4.5 mm
● NP (2.9 mmD)	3.7 mmD	CITNP31	CITNP33	CITNP34
● NP (2.9 mmD)	4.5 mmD	CITNP41	CITNP43	CITNP44
Replacement Retaining Screw		CASLT	CASLT	CASLT

Direct Transfers (Open-Tray Procedure) Color-coded by implant platform.

Includes a Retaining Screw (CASLC)



Catalog Numbers				
Implant Platform	Emergence Profile	Cuff Height		
		1.5 mm	3.0 mm	4.5 mm
● NP (2.9 mmD)	3.7 mmD	CDTNP31	CDTNP33	CDTNP34
● NP (2.9 mmD)	4.5 mmD	CDTNP41	CDTNP43	CDTNP44
Replacement Retaining Screw		CASLC	CASLC	CASLC

Implant Analog, Titanium Color-coded by implant platform.



Catalog Numbers	
Implant Platform	Catalog No.
● NP (2.9 mmD)	CIANP

Figure A



NP Platform

Figure B

The top surface of the healing collar is etched with three symbols to reference implant platform diameter (left), emergence profile diameter (top right) and cuff height (lower right). In some cases, only the initial digit of the measurement is shown. See item number table above for specific measurements. NP = Narrow Platform (2.9 mmD)



Titanium Temporary Abutment

Includes a Retaining Screw (CUAS) and a Long Processing Screw (CASLC)



Catalog Numbers		
Implant Platform	Emergence Profile	Catalog No.
● NP (2.9 mmD)	3.5 mmD	CTANP31
Replacement Retaining Screw		CUAS
Replacement Long Processing Screw		CASLC

Contour Abutments, Straight

Includes a Retaining Screw (CUAS)



Catalog Numbers				
Implant Platform	Emergence Profile	Cuff Height		
		1.5 mm	3.0 mm	4.5 mm
● NP (2.9 mmD)	3.7 mmD*	CANP31S	CANP33S	CANP34S
● NP (2.9 mmD)	● 4.5 mmD**	CANP41S	CANP43S	CANP44S
Replacement Retaining Screw		CUAS	CUAS	CUAS

* Contour Abutments with 3.7 mmD Emergence Profile are not compatible with 3.5 mmD Contour Restorative Components.

** Please refer to the Tapered Screw-Vent Implant System catalog for a list of 4.5 mmD Contour Restorative Components available to facilitate restoration of Contour Abutments with 4.5 mmD Emergence Profile.

Angled Abutments, 20°

Includes a Retaining Screw (CUAS)



Catalog Numbers		
Implant Platform	Emergence Profile	Catalog No.
● NP (2.9 mmD)	3.5 mmD	C20A3
Replacement Retaining Screw		CUAS

"Cast To" Gold Abutments, Engaging

Includes a Retaining screw (CUAS)



Catalog Numbers		
Implant Platform	Emergence Profile	Catalog No.
● NP (2.9 mmD)	3.5 mmD	CEANP31
Replacement Retaining Screw		CUAS
Long Processing Screw		CASLC

Contour Abutments, Angled 17°

Includes a Retaining Screw (CUASA or CUAS).



Catalog Numbers			
Implant Platform	Emergence Profile	Cuff Height	
		1.5 mm	3.0 mm
● NP (2.9 mmD)	3.7 mmD*	CANP31A	CANP33A
● NP (2.9 mmD)	4.5 mmD*	CANP41A	CANP43A
Replacement Retaining Screw		CUASA	CUAS

* Contour Abutments with 3.7 mmD Emergence Profile are not compatible with 3.5 mmD Contour Restorative Components.

** Please refer to the Tapered Screw-Vent Implant System catalog for a list of 4.5 mmD Contour Restorative Components available to facilitate restoration of Contour Abutments with 4.5 mmD Emergence Profile.

Locator® Abutments



Catalog Numbers						
Implant Platform	Cuff Height					
	1.0 mm	2.0 mm	3.0 mm	4.0 mm	5.0 mm	6.0 mm
● NP (2.9 mmD)	CLOCNP1	CLOCNP2	CLOCNP3	CLOCNP4	CLOCNP5	CLOCNP6

Locator Attachment Components Compatible with the existing Locator Components. Please refer to the Tapered Screw-Vent Catalog for additional options.



Catalog Numbers		
Description	Qty.	Catalog No.
Male Processing Package, Includes 2 Each of Metal Housings with Black Processing Males, Block-Out Spacers, Clear, Pink and Blue Replacement Males	2 Ea.	LOCMP2
Male Processing Package, Includes 10 Each of Parts Listed Above	10 Ea.	LOCMP10

Ball Abutment Abutments do not engage internal hex connection. Cap Attachment Housing and Nylon Liner are sold separately (CA).



Catalog Numbers				
Implant Platform		Cuff Height		
		2.0 mm	4.0 mm	6.0 mm
● NP (2.9 mmD)		CBANP2	CBANP4	CBANP6

Overdenture Components Compatible with the existing Overdenture Components. Please refer to the Tapered Screw-Vent Catalog for additional options.



Catalog Numbers		
	Description	Catalog No.
	Cap Attachment Housing (CAH)/Cap Attachment Nylon Liner (CAN)	CA
	Cap Attachment Housing	CAH
	Cap Attachment Nylon Liner (Transparent)	CAN
	Cap Attachment Nylon Liner (Gray — Rigid Retention)	CAN-G

Prosthetic Tools

Catalog Numbers		
	Description	Catalog No.
	Hex Driver, Short, with GemLock® Retention (1.25 mm, 22 mmL)	HXGR1.25
	Hex Driver, Long, with GemLock Retention (1.25 mm, 30 mmL)	HXLGR1.25
	Latch-Lock Hex Driver (Short) for Surgical and Prosthetic Screws (may require Reduction Handpiece) (1.25 mm, 23 mmL)	HX1.25D
	Latch-Lock Hex-Driver (Long) for Surgical and Prosthetic Screws (may require Reduction Handpiece) (1.25 mm, 26 mmL)	HXL1.25D
	Torque Wrench, Restorative (adjustable torque range 10 Ncm – 35 Ncm)	TWR
	Torque Wrench Hex Driver, Short (1.25 mm, 17 mmL)	TW1.25
	Torque Wrench Hex Driver, Long (1.25 mm, 22 mmL)	TW1.25L
	Removal Tool for Eztetic Abutments (NEW)	CLRT2


NOTE: CLRT2 is included with the NP Surgical Module.

Tapered Screw-Vent Surgical Kit




*Call a sales representative for availability in the kit.

NP Surgical Module



The NP Surgical Module snaps into the Tapered Screw-Vent Surgical Kit



For maximum cutting efficiency, replace drills frequently.

Color Coding

Implant diameter

3.1 mmD

Surgical sequence color bar



Drill band color for dense bone protocol



Implant cap color and restorative platform



Vial cap label

NP: Narrow Platform (2.9 mmD)




Can be inserted into
Tapered Screw-Vent
Surgical Kit (TSVKIT)

NP Surgical Module

Catalog Numbers		
Description	Qty.	Catalog No.
NP Surgical Module (To be inserted into the TSVKIT). Includes:	1 Ea.	NPMOD
NP Surgical Module (Tray Only)		NPTRAY
2.4 mmD Driva Surgical Drill, 22 mmL		SV2.4DN
2.8/2.4 mmD Driva Surgical Step Drill, 22 mmL		ZOP28DN
3.0 mmD Cortical Bone Tap Tool		ZOPTT30
2.1 mm GemLock Hex Drill, Latch Lock, Short		CHD2.1
2.1 mm GemLock Hex Drill, Latch Lock, Long		CHDL2.1
2.1 mm GemLock Hex Tool, Ratchet, Short		CHR2.1
2.1 mm GemLock Hex Tool, Ratchet, Long		CHRL2.1
Removal Tool for Eztetic Abutments		CLRT2

Please refer to the Tapered Screw-Vent Implant System Catalog for a complete list of surgical instrumentation in the Zimmer® Instrument Kit System.

Replacement Surgical Instrumentation




Catalog Numbers	
Description	Catalog No.
2.4 mmD Driva Surgical Drill, 22 mmL	SV2.4DN
2.8/2.4 mmD Driva Surgical Step Drill, 22 mmL	ZOP28DN
3.0 mmD Cortical Bone Tap Tool	ZOPTT30
2.1 mm GemLock Hex Drill, Latch Lock, Short	CHD2.1
2.1 mm GemLock Hex Drill, Latch Lock, Long	CHDL2.1
2.1 mm GemLock Hex Tool, Ratchet	CHR2.1
2.1 mm GemLock Hex Tool, Ratchet	CHRL2.1
Removal Tool for Eztetic Abutments	CLRT2

Drilling Sequence Guidelines

Soft bone protocol: follow solid color bars on the surgical tray surface until the segmented color bar. The segmented color bar indicates the final drill for soft bone protocol.


Dense bone protocol: follow solid color bars only. The last solid bar in the sequence represents the final drill for dense bone.



Surgical Procedure


Site Preparation

Please follow Site Preparation Instructions described in the Tapered Screw-Vent Implant System Surgical Manual along with a drilling sequence for the 3.1 mmD Eztetic Implant.




Removing The Implant From The Vial

Remove the implant outer vial from the box and open the outer vial to break the seal. Drop the sterile inner vial and contents onto a sterile field. Flip the white top of the inner vial open by pressing on the flat side with access hole. Press the top to the inner vial body to lock in the top. Grasp the top of the titanium packaging component placed on top of the implant, remove it and discard.



Delivering The Implant To The Site

Place the appropriate insertion instrument directly into the implant. The following instruments can be used for implant delivery to the site: the GemLock Hex Drill [CHD2.1, CHDL2.1] attached to a motor handpiece, or GemLock Hex Driver [CHR2.1, CHRL2.1] attached to the GemLock Retaining Square Ratchet [RSR] or Stainless Steel Screwdriver Handle [SSHS]. Carry the implant via the selected delivery instrument(s) to the receptor site and place directly into the prepared osteotomy.



Inserting And Orienting The Implant

Rotate the implant into place with the selected delivery instrument(s). The GemLock Hex Drills and Drivers are designed with six flats to align with the implant hex. To ensure proper orientation of the Contour Abutments, align the flat side of the Hex Drill or Driver to the buccal aspect. For Angled Abutments, orient a flat side of the Hex Drill or Driver toward the direction of the implant angle. Follow One-Stage or Two-Stage Healing Instructions provided in the Instructions for Use.

1. Data on file.
2. Chu C-M, Huang H-L, Hsu J-T, Fuh L-J. Influences of internal tapered abutment designs on bone stresses around a dental implant: three-dimensional finite element method with statistical evaluation. *J Periodontol* 2012;83:111-118.
3. Trisi P, Marcato C, Todisco M. Bone-to-implant apposition with machined and MTX microtextured implant surfaces in human sinus grafts. *Int J Periodontics Restorative Dent*. 2003;23(5):427-437.
4. Todisco M, Trisi P. Histomorphometric evaluation of six dental implant surfaces after early loading in augmented human sinuses. *J Oral Implantol*. 2006;32(4):153-166.
5. Shin SY, Han DH. Influence of a microgrooved collar design on soft and hard tissue healing of immediate implantation in fresh extraction sites in dogs. *Clin Oral Implants Res*. 2010;21:804-814.



This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

This image shows a full page of white paper with horizontal blue or grey ruling lines. The lines are evenly spaced and run across the width of the page, typical of notebook paper. There are no margins, text, or other markings on the page.



Contact us at 1-800-342-5454 or visit
zimmerbiometdental.com

Zimmer Biomet Dental
Global Headquarters
4555 Riverside Drive
Palm Beach Gardens, FL 33410
Tel: +1-561-776-6700
Fax: +1-561-776-1272

All references to Zimmer Biomet Dental contained herein refer to the dental subsidiaries of Zimmer Biomet Holdings, Inc. Unless otherwise indicated, as referenced herein, all trademarks are the property of Zimmer Biomet; and all products are manufactured, distributed and marketed by Zimmer Biomet Dental (and, in the case of distribution and marketing, its authorized marketing partners). Locator is a trademark of Zest IP Holdings LLC. Locator Attachments are manufactured by Zest Anchors, LLC. For additional product information, please refer to the individual product labeling/IFU. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. This material may not be copied or reprinted without the express written consent of Zimmer Biomet.

ZBINST1639 REV A 11/16 ©2016 Zimmer Biomet, All rights reserved.

