LMDental

LM-Instruments Oy PL 88 (Norrbyn rantatie 8) FI-21601 Parainen Finland

Tel. +358 2 4546 400 Fax +358 2 4546 444 info@lm-dental.com www.lm-dental.com

Contact your local dealer



Diamond like feel & precision - LM-ErgoSense Dark DIAMOND

LM-ErgoSense restorative instruments with Dark Diamond non-stick coating are ideal for placing and contouring of composite resins. The diamond like coating reduces the stickiness and allow for pleasant and precise restorative instrumentation.

Features:

- Exclusive diamond like coating
- Hard surface
- Smooth surface
- Dark color coating
- LM-ErgoSense handle

Benefits:

- Non-stick placement and sculpting of the composite resin
- Not sensitive to scratches
- Easy to clean
- Maximum contrast against the filling & minimal reflections
- Extreme precision and feel



Selection for maximum precision:

Burnisher LM 31-32DD ES

For packing and modeling composite.

Burnisher-Plugger, ball-ended 1.5-2.5 mm LM 471-473DD ES

For modeling composite.

Modeller MH LM 708-709DD ES

For packing large approximal fillings. Cylinder shaped working end will not get stuck into the composite while light curing.

Modeller MH Mini LM 707-708DD ES

For packing small and medium approximal fillings. Cylinder shaped working end will not get stuck into the composite while light curing.

Packer-Modeller IO LM 482-483DD ES

For packing and modelling composite. Conical working ends optimal for shaping occlusal surfaces.

Plastic instrument, anterior LM 43-47DD ES

For placing, packing and contouring composite.

Plastic instrument, posterior LM 48-702DD ES

For placing, packing and

contouring composite.

Mini Spatula LM 444-445DD ES

For placing and contouring composite.

Dark Diamond instruments can be cleaned and sterilized like any other LM instrument.