




Zfx™ GenTek™

Genuine Restorations and  
Open Digital Workflows

Product Catalog





# GenTek™ Open Connectivity

For even more flexibility.




## Genuine

GenTek™ Pre-Abutment Blanks are compatible with Medentika® holders. (CAM update required!)



## Compatible

Compatibility of the GenTek™ Ti-Base to Sirona® CEREC® Blocks



exocad

dental wings

3shape



## Free

Libraries for the GenTek restorative components are available for Zfx™ CAD, exocad, 3Shape and Dentalwings CAD software. Free download at [zfx-dental.com](http://zfx-dental.com)

#### Coming soon in 2019:

- Ti-Bases and ancillary components for Certain®, External Hex® and Eztec® connections
- Ti-Bases for full-arch restorations
- Pre-milled abutment blanks for all 4 connections



SCAN & INFORM  
More at  
[www.zfx-dental.com](http://www.zfx-dental.com)



A COMPANY OF

ZIMMER BIOMET

# Zfx™ GenTek™ Restorative Components

## Table of Contents

<b>Zimmer Biomet Dental Tapered Screw-Vent® and Trabecular metal connection</b>	<b>4</b>
<b>Zimmer Biomet Dental Eztetic® Connection</b>	<b>6</b>

## Placement Tool Compatibility Guide







Placement Tool	Laboratory Analog
M 16	Certain & Eztetic Lab Analogs
M 18	TSV, TM & Tapered Abutment Lab Analogs

## Legend

### Ø Restorative Platform (mm)

Color code for Eztetic® and Tapered Screw-Vent®

	2.9mmD
	3.5mmD
	4.5mmD
	5.7mmD

# Zfx™ GenTek™ Restorative Components: Implant-Level Connection

## Zimmer Biomet Dental Tapered Screw-Vent® and Trabecular Metal™ Connection

### Single-Unit Restorative Components

#### **Ti-Base Engaging<sup>1</sup>** incl. Retaining and Gluing Screws



Height (mm)	Restorative Platform (mm)	Part No.
4.7	3.5	ZFX11-ZB-TSV-3547-EL
	4.5	ZFX11-ZB-TSV-4547-EL
	5.7	ZFX11-ZB-TSV-5747-EL

#### **Retaining Screw**



Part No.	
ZFX09-ZB-TSV-RS	30 Ncm

#### **Titanium Pre-milled Abutment Blank<sup>2</sup>** incl. Retaining Screw



Restorative Platform (mm)	Part No.
3.5	ZFX18-ZB-TSV-35
4.5	ZFX18-ZB-TSV-45
5.7	ZFX18-ZB-TSV-57

#### **Retaining Screw**



Part No.	
ZFX09-ZB-TSV-RS	30 Ncm

### Multi-Unit Restorative Components

#### **Ti-Base Non-Engaging<sup>1</sup>** incl. Retaining and Gluing Screws



Height (mm)	Restorative Platform (mm)	Part No.
4.7	3.5	ZFX11-ZB-TSV-3547-NEL
	4.5	ZFX11-ZB-TSV-4547-NEL
	5.7	ZFX11-ZB-TSV-5747-NEL

#### **Retaining Screw**




Part No.	
ZFX09-ZB-TSV-RS	30 Ncm



## Components Accessories & Tools

### Desktop Scanbody incl. Screw



Restorative Platform (mm)	Part No.	Hand-tighten <sup>3</sup>
3.5	ZFX05D-ZB-TSV-35	
4.5	ZFX05D-ZB-TSV-45	
5.7	ZFX05D-ZB-TSV-57	

### Torque Wrench for Desktop Scanbody



Part No.
ZFX02002067

### Digital Scan Analog<sup>4</sup>

### Placement Tool for Digital Scan Analogs



Restorative Platform (mm)	Part No.	Part No.
3.5	ZFX15-ZB-TSV-35	ZFX02-M18
4.5	ZFX15-ZB-TSV-45	
5.7	ZFX15-ZB-TSV-57	

### Intraoral Scanbodies incl. Screws



Height (mm)	Restorative Platform (mm)	Part No.	Hand-tighten <sup>3</sup>
4.0 & 7.0	3.5	ZFX05I-ZB-TSV-35	
	4.5	ZFX05I-ZB-TSV-45	
	5.7	ZFX05I-ZB-TSV-57	

### Latch-Lock Hex Driver<sup>5</sup> for Abutment Screws



Description	Part No.
Latch-Lock Hex Driver for Surgical and Prosthetic Screws, Short, 1.25 mm D, 23 mm L (May require Reduction Handpiece)	HX1.25D
Latch-Lock Hex Driver for Surgical and Prosthetic Screws, Long, 1.25 mm D, 26 mm L (May require Reduction Handpiece)	HXL1.25D

### Hex Driver for Abutment Screws & Intraoral Scanbodies



Description	Part No.
Hex Driver, Short, with GemLock(R) Retention, 1.25 mm D, 22 mm L	HXGR1.25
Hex Driver, Long, with GemLock(R) Retention, 1.25 mm D, 30 mm L	HXLGR1.25

### Removal Tool



Description	Part No.
Removal Tool for TVS/TM Abutments	TLRT2
Removal Screw for TSV/TM Abutments	HLRS

### Torque Wrench Tool for Abutment Screws & Intraoral Scanbodies



Description	Part No.
Torque Wrench, Restorative (adjustable torque range 10 – 35 Ncm)	TWR
Torque Wrench Driver, Short	TW1.25 <sup>5</sup>
Torque Wrench Driver, Long	TW1.25L <sup>5 6</sup>

<sup>1</sup>Groove indicating 4 mm post height

<sup>2</sup>Coming soon

<sup>3</sup>Minimum of 5Ncm; Maximum of 15Ncm.  
Refer to Instruction Manual for more information.

<sup>4</sup>Compatible with existing TSV/ TM Impression Copings

<sup>5</sup>Compatible with Tibases

<sup>6</sup>Compatible with IOS Scanbodies

# Zfx™ GenTek™ Restorative Components: Implant-Level Connection

## Zimmer Biomet Dental Eztetic® Connection

### Single-Unit Restorative Components

#### **Ti-Base Engaging<sup>1</sup>** *incl. Retaining and Gluing Screws*



Height (mm)	Restorative Platform (mm)	Part No.
4.7	2.9	ZFX11-ZB-EZ-2947-ES

#### **Retaining Screw**



Part No.	Recommended Torque
ZFX09-ZB-EZ-RS	30 Ncm

#### **Titanium Pre-milled Abutment Blank** *incl. Retaining Screw*



Restorative Platform (mm)	Part No.
2.9	ZFX18-ZB-EZ-29

#### **Retaining Screw**



Part No.	Recommended Torque
ZFX09-ZB-EZ-RS	30 Ncm

### Multi-Unit Restorative Components

#### **Ti-Base Non-Engaging<sup>1</sup>** *incl. Retaining and Gluing Screws*



Height (mm)	Restorative Platform (mm)	Part No.
4.7	2.9	ZFX11-ZB-EZ-2947-NES




#### **Retaining Screw**







Part No.	Recommended Torque
ZFX09-ZB-EZ-RS	30 Ncm



## Components Accessories & Tools

Desktop Scanbody <i>incl. Screw</i>			Torque Wrench <i>for Desktop Scanbody</i>		Digital Scan Analog <sup>3</sup>		Placement Tool <i>for Digital Scan Analogs</i>	
								
Restorative Platform (mm)	Part No.		Part No.		Restorative Platform (mm)	Part No.		Part No.
2.9	ZFX05D-ZB-EZ-29	Hand-tighten <sup>2</sup>	ZFX02002067		2.9	ZFX15-ZB-EZ-29	ZFX02-M16	

Intraoral Scanbodies <i>incl. Screws</i>			Latch-Lock Hex Driver <sup>4</sup> <i>for Abutment Screws</i>	
				
Height (mm)	Restorative Platform (mm)	Part No.	Description	Part No.
4.0 & 7.0	2.9	ZFX05I-ZB-EZ-29	Latch-Lock Hex Driver for Surgical and Prosthetic Screws, Short, 1.25 mm D, 23 mm L May require Reduction Handpiece)	HX1.25D
			Latch-Lock Hex Driver for Surgical and Prosthetic Screws, Long, 1.25 mm D, 26 mm L (May require Reduction Handpiece)	HXL1.25D

Hex Driver <i>for Abutment Screws &amp; Intraoral Scanbodies</i>		Removal Tool		Torque Wrench Tool <i>for Abutment Screws &amp; Intraoral Scanbodies</i>	
					
Description	Part No.	Description	Part No.	Description	Part No.
Hex Driver, Short, with GemLock(R) Retention, 1.25 mm D, 22 mm L	HXGR1.25	for Eztec Abutments	CLRT2	Torque Wrench, Restorative (adjustable torque range 10–35 Ncm)	TWR
Hex Driver, Long, with GemLock(R) Retention, 1.25 mm D, 30 mm L	HXLGR1.25			Torque Wrench Driver, Short	TW1.25 <sup>4</sup>
				Torque Wrench Driver, Long	TW1.25L <sup>4 5</sup>

<sup>1</sup>Groove indicating 4 mm post height

<sup>2</sup>Eztetic connection TI-Base require the use of the CEREC blocks for Certain 3.4 mm Connection

<sup>3</sup>Minimum of 5Ncm; Maximum of 15Ncm. Refer to Instruction Manual for more information.

<sup>4</sup>Compatible with existing Eztec Impression Copings

<sup>5</sup>Compatible with Tibases

<sup>6</sup>Compatible with IOS Scanbodies



To order genuine  
GenTek™ Components,  
call: +34 93 470 55 00

**Zfx GmbH**

Kopernikusstraße 15  
85221 Dachau, Germany  
Tel. +49 81 31 33 244 - 0  
Fax +49 81 31 33 244 - 10  
[www.zfx-dental.com](http://www.zfx-dental.com)

